



Introducing...
THE BOOST 2
Effortless Unweighting



Combining Boost's trusted unweighting functionality with new enhancements designed to make usability & serviceability effortless.

BOOST 2 FEATURES



The Most Sophisticated Micro-Gravity Treadmill™ on the Planet.



PRECISION HEIGHT ADJUSTMENT

Our motorized height adjustment lifts precisely from 4'6" to 7'. With no manual lifting involved, it makes for safe, quick and easy exit/entry.



SIMPLIFIED USER CONTROLS

No more pressing and holding buttons! Our new control system is easy and intuitive, featuring knob controls that allow for quick and accurate adjusting of speed and body weight.



SPEED

Up to 18mph



CLINICAL GAIT VIEWS

Front, side & rear



INCLINE

Up to 15%



WHISPER QUIET

Smooth & silent



FOOTPRINT

82" L x 42" W



TREADMILL

Woodway® 4FRONT



SHORTS

Thin, comfortable,
& durable



WEIGHT-BEARING

20-100% body weight
adjustment